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Manufacturer and US Distributor: Fresenius Medical Care North America (FMCNA), 920 Winter St. Waltham, MA 02451, 800-323-5188

Canadian Importer/Distributor: Fresenius Medical Care Canada Inc., 45 Staples Avenue, Suite 110, Richmond Hill, Ontario L4B 4W6 888-709-4411

Imported/Distributed in Mexico by: Fresenius Medical Care de México S.A. de C.V., Paseo del Norte 5300, "Guadalajara Technology Park", CP 45010, Zapopan, Jal. México (52) (33) 35 40 42 00

Printed in U.S.A. 71-4400.00 07/12

Instructions for Use
650g and 900g *bi*bag® Bicarbonate Concentrate

Powder for production of liquid sodium bicarbonate concentrate for bicarbonate hemodialysis

INDICATIONS FOR USE: The Fresenius Medical Care North America (FMCNA) *bi*bag system is used with FMCNA three stream proportioning hemodialysis systems equipped with the *bi*bag module, such as the 2008T Hemodialysis Machine, and is intended for use in bicarbonate hemodialysis for acute and chronic renal failure. *bi*bag bicarbonate concentrate is intended for extracorporeal bicarbonate hemodialysis according to a physician's prescription.

CONTRAINDICATIONS: There are no absolute contraindications for bicarbonate dialysis. *bi*bag bicarbonate concentrate must be used only on the direction of a physician who has evaluated all the pertinent features of the product in relation to the individual patient.

REQUIREMENTS:

- For use only with FMCNA dialysis machines with three-stream proportioning bicarbonate systems when set to proportion 1 part acid to 44 parts bicarbonate base.
- For use only with FMCNA 45x acid concentrates.
- For use only with water that meets current ANSI/AAMI RD62, Ph Eur or ISO 13959.
- For use only with FMCNA systems equipped with the *bi*bag connector and Diasafe Plus filter.

WARNING:

- Failure to follow these Instructions for Use may cause patient injury or death.
- NOT FOR PARENTERAL USE. For hemodialysis use only.
- Do not use damaged bag.
- Product must be used before the end of expiration month (Exp. Date) on the product's label.
- The dissolved *bi*bag bicarbonate concentrate must be used within 24 hours of connecting to the dialysis machine.
- Always check the composition of the dialysate solution before use.
- Dispose of the bag and unused powder according to standard unit protocols.
- Dispose of the bag and unused powder between each treatment.
- Use in accordance with machine operating instructions.
- Bacterial growth may occur when using bicarbonate concentrate.

CAUTION:

- Federal (USA) law restricts this product to sale by, or on order of, a physician.
- Prevent any damage to the filter at the end of the tube inside the bag. Do not use *bi*bag bicarbonate concentrate if this filter is broken.
- The information in this Instructions for Use must be carefully read and understood before using the product.
- Do not refill the *bi*bag container.
- The operation of the specially adapted dialysis machine is restricted to persons who are trained in dialysis and who are completely familiar with the operating instructions stated in the operator's manual of the corresponding dialysis machine and with these Instructions for Use.

ADVERSE REACTIONS: Some undesirable effects related to dialysis can occur, for instance, hypotension, nausea, headache and cramps.

DIRECTIONS FOR USE:

1. Remove seal from *bi*bag bicarbonate concentrate and inspect the product for visible damage. If the bag or bag top is damaged replace with a new *bi*bag.
2. Open the *bi*bag door on the machine by lifting up on the dark gray handle.
3. With the white handle facing outward, hang the *bi*bag on the *bi*bag connector nozzles.
4. Close door of connector, making sure it latches firmly in place. An audible click will indicate the door is closed.
5. At the end of treatment, the *bi*bag may be emptied, for disposal, either manually or by using the "Empty *bi*bag" function of the dialysis machine.
6. Refer to the detailed directions for use provided in the dialysis machine operator's manual.

RUN TIME: The *bi*bag disposable contains a fixed volume of bicarbonate powder. Refer to Table 1 below to verify that enough run time is available to complete a treatment using one bag. *Run times are estimates and may vary per unit protocol. Run times include estimated preparation and treatment time. Highlighted run times indicate 10+ hours of treatment time. If desired bicarbonate or flow setting not found, please refer to the Hemodialysis Machine *bi*bag System Operator's Instructions.

Table 1: <i>bi</i> bag Disposable Estimated Run Time (Minutes)*						
		Bicarbonate Setting (mEq/L or mmol/L)				
bi		40	36	32	28	
Dialysate flow QD (mL/min)	650g	800	180	200	225	257
		700	206	229	257	294
		600	240	267	300	343
		500	288	320	360	411
		400	360	400	450	514
		300	480	533	600	686
	900g	200	720	800	900	1029
		800	240	267	300	343
		700	274	305	343	392
		600	320	356	400	457
		500	384	427	480	549
		400	480	533	600	686
		300	640	711	800	914
		200	960	1067	1200	1371

STORAGE: Store between +15°C to +30°C. Do not refrigerate. In order to avoid any contamination, keep the bag tightly sealed until it is connected to the hemodialysis machine.

CHEMICAL COMPOSITION OF CONTENT: Sodium Bicarbonate, USP/Ph Eur.

PRESENTATION: 650g/900g *bi*bag bicarbonate concentrate in polyamide/polyethylene container.

NOTE: In some cases the bicarbonate powder may contain particles that have been slightly discolored after exposure to water. These discolorations originate from the raw material (sodium bicarbonate). Both the powder and the liquid solution obtained from these *bi*bag products meet the demands of the European and US Pharmacopoeia. The safety and efficacy of the *bi*bag product are not affected by these discolorations.

Printed In the USA 71-4400.00 07/12

Combined
*bi*bag Instructions (Inside)

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